

**Anti-PIVKA-II antibody (23-46aa) [12103] (STJ400349)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	CLIA/FIA
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

**TARGET INFORMATION**

Immunogen Region	23-46aa
Specificity	Antibody recognizes human PIVKA-II (Protein Induced Vitamin K absence or Antagonist-II). Recognizes Prothrombin.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---