

Anti-Cannabinoid antibody (STJ400355)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	THC

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	≥ 1.5 mg/ml
Conjugation	Unconjugated
Formulation	0.05 M Potassium Phosphate Buffer, pH 7.2 with 0.85 % NaCl, 0.05 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes THC and its metabolites. Does not react with Opiates, Benzoylcegonine, Amphetamines, Phencyclidine and Cannabidiol.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---