

Anti-Clonazepam antibody (STJ400368)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Clonazepam

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	≥ 1.5 mg/ml
Conjugation	Unconjugated
Formulation	0.05 M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody primarily recognizes clonazepam. 7-aminoclonazepam, 7-aminoflunitrazepam, Delorazepam >10%; N-Desalkylflurazepam, Flunitrazepam >1%; Nitrazepam, Desmethyldiazepam, Diazepam
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---