

## Anti-hCG beta antibody [5011] (STJ400380)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CLIA/FIA
Host / Source	Mouse
Reactivity	hCG

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	29 mM citrate, 142 mM phosphate, pH 6.5, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Antibody recognizes human chorionic gonadotropin and its free beta subunit. Does not recognize LH
-------------	---

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---