

## Anti-GFAP antibody (60-383aa) [R13501] (STJ400402)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FIA                |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 5.0 mg/ml (+/-10 %)   |
| Conjugation         | Unconjugated  |
| Formulation         | 50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative |
| Isotype             | IgG1  |
| Storage Instruction | Store at 2–8 °C, and do not freeze.   |

### TARGET INFORMATION

|                  |   |
|------------------|---|
| Immunogen Region | 60-383aa  |
| Specificity      | Antibody recognizes human GFAP. Does not recognize vimentin, desmin or peripherin |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|