

Anti-GFAP antibody (60-383aa) [R13503] (STJ400404)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FIA |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 5.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Formulation | 50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| Isotype | IgG1 |
| Storage Instruction | Store at 2–8 °C, and do not freeze. |

TARGET INFORMATION

| | |
|------------------|---|
| Immunogen Region | 60-383aa |
| Specificity | Antibody recognizes human GFAP. Does not recognize desmin or peripherin. Possibility to a minor cross-reaction to vimentin. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|