

Anti-GFAP antibody (60-383aa) [R13504] (STJ400405)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Immunogen Region	60-383aa
Specificity	Antibody recognizes human GFAP. Does not recognize desmin or peripherin. Possibility to a minor cross-reaction to vimentin.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---