

## Anti-Amyloid beta antibody (34-40aa) [R13611] (STJ400411)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FIA
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	50 mM citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Immunogen Region	34-40aa
Specificity	Antibody recognizes the C-terminus of amyloid beta peptide 1-40. Does not recognize amyloid beta 1-38 or 1-42

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------