

Anti-HIV-1 p24 antibody [R14101] (STJ400414)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/LF/ELISA/FIA
Host / Source	Rabbit
Reactivity	HIV-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline (PBS) pH 7.4 with 0.095 % NaN ₃ as a preservative
Isotype	IgG
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Antibody recognizes human immunodeficiency virus 1 p24 (HIV-1 p24). Recognizes HIV-1 p24 group O and N. Antibody also recognizes HIV-1 p24 group M clades A, B, C and D.
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---