

Anti-CD18 antibody [MEM48] (STJ400438)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/FC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Recognizes the CD18 molecule, a 95 kDa integrin molecule ($\beta 2$ integrin) that forms the common beta chain of the CD11/CD18 heterodimer integrins..
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---