

Anti-CD29 antibody [B-D15] (STJ400439)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC/FC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Formulation | Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. |
| Isotype | IgG2ak |
| Storage Instruction | Store at 2–8 °C, and do not freeze. |

TARGET INFORMATION

| | |
|-------------|--|
| Specificity | Recognises the GPIIa, integrin b1, VLA-b chain, a 130 kDa protein. |
|-------------|--|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|