

## Anti-CD178 antibody [B-R17] (STJ400467)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/FC/FUNC      |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1.0 mg/ml (+/-10 %)  |
| Conjugation         | Unconjugated   |
| Formulation         | Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. |
| Isotype             | IgG1k  |
| Storage Instruction | Store at 2–8 °C, and do not freeze.  |

### TARGET INFORMATION

|             |  |
|-------------|--|
| Specificity | Recognises APO1L, Fas ligand (FAS L) , a 40 kDa protein. |
|-------------|--|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|