

## Anti-CD31 antibody [B-N14] (STJ400485)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Isotype	IgG1k
Storage Instruction	Store at 2-8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Recognises the Platelet Endothelial Cell Adhesion Molecule-1 (PECAM-1) , GPIIa, a 130-140 kDa protein.
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---