

Anti-CD141 antibody [B-S24] (STJ400499)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 1.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Formulation | Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. |
| Isotype | IgG1k |
| Storage Instruction | Store at 2–8 °C, and do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------------------------------|
| Specificity | Recognises thrombomodulin (TM) , a 105 kDa protein. |
|-------------|-----------------------------------------------------|

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|