

Anti-CD261 antibody [B-T35] (STJ400503)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Formulation | Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. |
| Isotype | IgG1k |
| Storage Instruction | Store at 2–8 °C, and do not freeze. |

TARGET INFORMATION

| | |
|-------------|---|
| Specificity | Recognises both natural and recombinant TRAIL-R1 (DR4) , a 43-46 kDa protein. |
|-------------|---|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|