

## Anti-CD1b/c antibody [B-B5] (STJ400513)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                              |
|---------------------|--------------------------------------------------------------|
| Clonality           | Monoclonal                                                   |
| Conjugation         | Unconjugated                                                 |
| Formulation         | Phosphate-buffered saline with 1% BSA and 0.09% sodium azide |
| Isotype             | IgG1k                                                        |
| Storage Instruction | Store at 2–8 °C, and do not freeze.                          |

### TARGET INFORMATION

|             |                                 |
|-------------|---------------------------------|
| Specificity | Recognises a 43-45 kDa protein. |
|-------------|---------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|