

Anti-CD21 antibody [B-E5] (STJ400537)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/FUNC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Isotype	IgG2ak
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Recognises the Epstein-Barr virus receptor (EBV R) , CR2, C3d R, a 110 kDa protein.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---