

Anti-CD59 antibody [MEM43] (STJ400556)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/FC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Isotype	IgG2ak
Storage Instruction	Store at 2-8 °C, and do not freeze.

TARGET INFORMATION

Specificity	This antibody reacts with a PI-linked glycoprotein, molecular weight NA 18-20 kDa, found in all types of leucocytes including platelets..
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---