

## Anti-CD218a antibody [B-E43] (STJ400592)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Formulation         | Phosphate-buffered saline with 1% BSA and 0.09% sodium azide |
| Isotype             | IgG1k  |
| Storage Instruction | Store at 2–8 °C, and do not freeze.                          |

### TARGET INFORMATION

|             |  |
|-------------|--|
| Specificity | Recognises the IL-18Ra, a 63 kDa protein. Cross-reacts with mouse CDw218a protein in cytometry.. |
|-------------|--|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|