

## Anti-CD8 antibody [B-Z31] (STJ400597)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Recognises the T-cell coreceptor, a 32-34 kDa protein.
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---