

Anti-CD112 antibody [B-C12] (STJ400674)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Isotype	IgG2bk
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Recognises Nectin-2, the polio retrovirus receptor-2 (PRR2).
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---