

Anti-HLA Class I antibody [W6/32] (STJ400706)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Recognises an antigenic determinant shared among products of the HLA-A,-B and-C loci.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---