

## Human IL-10 (IL-10) Antibody Pair (STJ400798)

### GENERAL INFORMATION

|               |                |
|---------------|----------------|
| Product Type  | Antibody Pairs |
| Applications  | ELISA/ELISpot  |
| Host / Source | Mouse          |
| Reactivity    | Human          |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1.0 mg/ml (+/-10 %) /0.1 mg/ml (+/-10 %)   |
| Conjugation         | Unconjugated/Biotin  |
| Formulation         | Capture: Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.<br>Detection: Phosphate-buffered saline with 1% BSA and 0.09% sodium azide |
| Isotype             | IgG1k  |
| Storage Instruction | Store at 2–8 °C, and do not freeze.  |

### TARGET INFORMATION

|             |  |
|-------------|--|
| Specificity | Recognises both natural and recombinant human IL-10./Recognises both natural and recombinant human IL-10 |
|-------------|--|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|