

Human CD62E (CD62E) Antibody Pair (STJ400803)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Concentration	1.0 mg/ml (+/-10 %) /0.1 mg/ml (+/-10 %)
Conjugation	Unconjugated/Biotin
Formulation	Capture: Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. Detection: Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Isotype	IgG1L/IgG1k
Storage Instruction	Store at 2–8 °C, and do not freeze./Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Recognises the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1) , E-selectin, a 97 kDa protein./Recognises the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1) , E-selectin, a 97 kDa protein.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS./STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---