

Human p53 (p53) Antibody Pair (STJ400804)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Concentration	1.0 mg/ml (+/-10 %) /0.1 mg/ml (+/-10 %)
Conjugation	Unconjugated/Biotin
Formulation	Capture: Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. Detection: Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Isotype	IgG2bk/IgG2ak
Storage Instruction	Store at 2–8 °C, and do not freeze./Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Recognises wild type and mutant p53, the epitope is within residues 18-30 of human p53./Recognises wild type and mutant p53, the epitope is within residues 18-30 of human p53.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS./STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---