

## Human IL-17F (IL-17F) Antibody Pair (STJ400808)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA/ELISpot
Host / Source	Mouse/Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Concentration	1.0 mg/ml (+/-10 %) /0.1 mg/ml (+/-10 %)
Conjugation	Unconjugated/Biotin
Formulation	Capture: Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free. Detection: Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Isotype	IgG2bk/IgG2bk
Storage Instruction	Store at 2–8 °C, and do not freeze./Store at 2–8 °C, and do not freeze.

### TARGET INFORMATION

Specificity	Recognises both natural and recombinant IL-17F./Recognises both natural and recombinant IL-17F.
-------------	---

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS./STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---