

Anti-Goat anti-Mouse Immunoglobulin G antibody (STJ400844)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Goat
Reactivity	IgG

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity Purified Polyclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.2
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	All IgG subclasses and all mouse Ig.
-------------	--------------------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---