

Anti-Follicle Stimulating Hormone Beta antibody (STJ400851)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	FSH Beta

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	0.05 M Phosphate Buffered Saline, pH 7.4
Isotype	IgG
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	FSH. FSH 100%, β hCG
-------------	----------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---