

Anti-Adenovirus antibody (STJ400858)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/LF/ELISA/IFA
Host / Source	Mouse
Reactivity	Adenovirus

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.0-7.5
Isotype	IgG2a
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Hexon Protein.
-------------	----------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---