

Anti-Herpes Simplex Virus 1 antibody (STJ400866)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IFA
Host / Source	Mouse
Reactivity	HSV 1-gE

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG2ak
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	gE.
-------------	-----

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---