

## Anti-Herpes Simplex Virus 2 antibody (STJ400867)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA/IFA       |
| Host / Source | Mouse              |
| Reactivity    | HSV 2-gG           |

### PRODUCT PROPERTIES

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                   |
| Conjugation         | Unconjugated                                                                                 |
| Purification        | Purified Monoclonal Antibody                                                                 |
| Formulation         | 10 mM Phosphate Buffered Saline, pH 7.4                                                      |
| Isotype             | IgG1k                                                                                        |
| Storage Instruction | Short Term ( $\leq 2$ weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |     |
|-------------|-----|
| Specificity | gG. |
|-------------|-----|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|