

Anti-Fentanyl antibody (STJ400872)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF
Host / Source	Mouse
Reactivity	Fentanyl

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified Monoclonal Antibody
Formulation	Phosphate Buffered Saline, pH 7.4
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Fentanyl. Norfentanyl 5%, Carfentanil oxalate 0.1%
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---