

Anti-African Swine Fever Virus antibody (STJ400928)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	ASFV-pp62

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Ion Exchange Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Isotype	IgG2b
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	ASFV pp62 protein.
-------------	--------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---