

Anti-Cancer Antigen 15-3 antibody (STJ400944)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CLIA
Host / Source	CHO Cells
Reactivity	CA 15-3

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffer, pH 7.4 with 0.15M NaCl
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---