

Anti-Influenza A Virus antibody (STJ400966)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF
Host / Source	CHO Cells
Reactivity	Influenza A Virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG2a
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Influenza A Virus Nucleoprotein. No reactivity with Influenza B, Adenovirus, RSV, Parainfluenza, M. pneumoniae or SARS-CoV-2
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---