

## Anti-Benzoylecgonine antibody (STJ400992)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	BE

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.2
Isotype	IgG
Storage Instruction	Short Term ( $\leq$ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	Benzoylecgonine and Cocaine.
-------------	------------------------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---