

Anti-Aripiprazole antibody (STJ401034)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Aripiprazole

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate Buffer with 0.85% NaCl, pH 7.2
Isotype	IgG1
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Aripiprazole. Aripiprazole 100%; Clozapine, Olanzapine, Quetiapine, Risperidone, Ziprasidone, Amitriptyline, Clomipramine, Cyclobenzaprine, Desipramine, Doxepin, Imipramine, Maprotiline, Nordoxepin, Nortriptyline, Trimipramine
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---