

Anti-Carbapenemase antibody (STJ401091)

GENERAL INFORMATION

Product Type	Primary antibodies
Host / Source	Mouse
Reactivity	NDM

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified Monoclonal Antibody
Formulation	Phosphate Buffered Saline, pH 7.4
Isotype	IgG1
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---