

Anti-Distemper Virus antibody (STJ401099)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline
Isotype	IgG2a
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	CDV surface envelope antigen.
-------------	-------------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---