

## Anti-Hepatitis B Surface Antigen antibody (STJ401119)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	HBsAg

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Ion Exchange Purified Monoclonal Antibody
Formulation	20 mM PBS pH 7.4, 50 mM NaCl, 0.1% NaN <sub>3</sub>
Isotype	IgG1
Storage Instruction	Short Term ( $\leq$ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	HBsAg (including adw, ayw, and adr subtypes).
-------------	---

### ADDITIONAL INFORMATION

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.</b>
------	--