

Anti-Ferritin antibody (STJ401136)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Ferritin

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	DEAE Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Isotype	IgG2
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Human Ferritin.
-------------	-----------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---