

## Anti-Benzodiazepine antibody (STJ401174)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Benzodiazepine (Oxazepam)

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Isotype	IgG1k
Storage Instruction	Short Term ( $\leq$ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	Oxazepam and related compounds including Alprazolam, Diazepam, Lorazepam, Nitrazepam and Tremazepam. Minimal reactivity with Clonazepam, Flurazepam, Prazepam
-------------	---

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---