

## Anti-Clostridium difficile toxin B antibody (STJ401204)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Clostridium difficile toxin B

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	Phosphate Buffered Saline, pH 7.0
Isotype	IgG1
Storage Instruction	Short Term ( $\leq$ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	C. difficile toxin B.
-------------	-----------------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---