

## Anti-Caffeine antibody (STJ401234)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | LF/ELISA           |
| Host / Source | Mouse              |
| Reactivity    | Caffeine           |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Purification        | Ion Exchange Purified Monoclonal Antibody  |
| Formulation         | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl                                    |
| Isotype             | IgG1k  |
| Storage Instruction | Short Term ( $\leq 2$ weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |           |
|-------------|-----------|
| Specificity | Caffeine. |
|-------------|-----------|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|