

Anti-Dengue Virus antibody (STJ401263)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF
Host / Source	Mouse
Reactivity	DENV-pan NS1

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffer, pH 7.4 with 0.15M NaCl
Isotype	IgG1
Storage Instruction	Store at-20 °C and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	NS1 protein of Dengue Virus types 1, 2, 3, 4.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---