

Anti-Hemoglobin A1c antibody (STJ401296)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/CLIA/ELISA
Host / Source	Mouse
Reactivity	HbA1c

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG1
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---