

Anti-Chikungunya antibody (STJ401307)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	CHIKV-Env 2

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Isotype	IgG2ak
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	CHIKV E2. No reactivity with Dengue, Zika, YFV, WNV, TBV, JEV
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---