

Anti-Helicobacter pylori antibody (STJ401330)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA/TIA
Host / Source	Mouse
Reactivity	H.pylori

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG1
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	H. pylori. No reactivity with Campylobacter, Clostridium difficile, E. coli, L. monocytogenes, Salmonella, Shigella, Y. enterocolitica
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---