

Anti-Adenovirus Type 1 & 2 antibody (STJ401338)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ELISA/IFA
Host / Source	Mouse
Reactivity	Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.2
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Canine Adenovirus Types 1 & 2. No reactivity with Canine Coronavirus, Canine Distemper Virus, Canine Parvovirus, Canine Parainfluenza Virus
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---