

Anti-Aldosterone antibody (STJ401344)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Aldosterone

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Ion Exchange Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Isotype	IgG2b
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Aldosterone. 17a-Hydroxyprogesterone, Estrone, 17a-Ethynylestradiol, 17 β -Estradiol, Estriol, Testosterone, Testosterone glucuronide, 5a-Dihydrotestosterone, Epitestosterone, 6 β -Hydroxytestosterone, 17a-Methyltestosterone, Progesterone, Androstenedion
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---